



Networked Robotics Corp.
825 Chicago Ave, Suite F
Evanston, IL 60202, USA
Toll free: (877) FRZ-TEMP
(877) 379-8367

info@networkedrobotics.com
support@networkedrobotics.com
NetworkedRobotics.com

Release Notes - Beta International Versions of the Tempurity Monitor for Windows 7,8,10

July 2017

Introduction

The following language versions of the Tempurity Monitor are available in this release. Because of limitations in alarm notification, these versions are of best use for watching the status of a population of “monitored devices”. All versions now run on Windows 7, 8, and 10.

Chinese	Hebrew*	Russian*
Danish*	Italian	Spanish
Dutch	Japanese	Swedish*
French	Malay*	Thai*
German	Portuguese	Vietnamese*

*=new in 2017

Testing

After downloading, enter a Tempurity Server name of “demo.networkedrobotics.com”.

Notes and Limitations

-In these beta releases of international language versions of the Tempurity Monitor alarm notifications are sent in English. Test alarm notifications are sent in the language specified.

-Voice phone alarm notifications are active in English only.

-These versions have not yet been validated for FDA GxP regulatory compliance.

-Some terms may remain in English.

-The names of “monitored devices” and some standard terms remain in the language used by the Tempurity Server to which the Tempurity Monitor is connected.

Change History

Previously several international versions were available for Windows XP platforms. This upgrade expands the number of languages supported by the Tempurity Monitor to 16 and makes all versions compatible with Windows 7, 8, 10.

Previous versions ran automatically in the spoken/written language that was set in Windows. These versions run in the language selected regardless of Windows setting.